

Computing on TSCC

- Make a folder for the class and move into it
 - `mkdir -p /oasis/tsc/scratch/username/biom262_harismendy`
 - `cd /oasis/tsc/scratch/username/biom262_harismendy`
- Make symbolic link to data folder
 - `ln -s /projects/ps-yeolab/biom262_2016/data/harismendy_data`
- Adjust your environment to have access to tools
 - `Source exportPATH.txt`
- Use screen
 - `screen -S name` to create
 - `^A ^D` to detach
 - `Screen -x [PID].name` to attach
 - `exit` to kill

UNIX commands

- Count number of lines
– `wc -l TSG.bed`
- Select lines with "TP53"
– `grep 'TP53' TSG.bed`
- Sort according to column 4
– `sort -k 4 TSG.bed`
- Select genes on chromosome 1
– `awk '$1=="chr1"' TSG.bed`
- Calculate gene length
– `awk '$5=$3-$2' TSG.bed`

Working with Intervals

Look at the file

```
more harismendy_data/class/CGC.exons.bed
```

Count number of lines

```
wc -l harismendy_data/class/CGC.exons.bed
```

Count the number of genes

```
cut -f 4 harismendy_data/class/CGC.exons.bed | sort | uniq  
| wc -l
```

Count number of intervals per gene

```
cut -f 4 harismendy_data/class/CGC.exons.bed | sort | uniq  
-c | sort -r | more
```

Calculate the size of each interval and sort by size

```
awk '{size=$3-$2; print $0,size}'  
harismendy_data/class/CGC.exons.bed | sort -n -r -k 5
```

Calculate the average interval size

```
awk '{sum=sum+$3-$2} END {print sum/NR}'  
harismendy_data/class/CGC.exons.bed
```

fastqc

Inspect the fastq file

```
zcat harismendy_data/class/file.fastq.gz | more
```

Read the help menu

```
fastqc -h
```

Create a output directory

```
mkdir fastqc
```

Run fastqc

```
fastqc -o fastqc  
harismendy_data/class/CGC.exons.bed/file.fastq.gz
```

Transfer the results to your desktop

Open results with web-browser

BWA alignment and BAM files

Read the doc

bwa mem or <https://github.com/lh3/bwa>

Start a Screen

screen -S username

Alignment + convert to sorted bam

bwa mem harismendy_data/resources/hg19_lite.fasta 1 R1.fq.gz R2.fq.gz | samtools view -bSh - > sample.bam

Sort and index the bam file

**samtools sort -m 2G sample.bam sample.sorted
samtools index sample.sorted.bam**

Remove Duplicate reads

- `java -jar`
`/opt/biotoools/picard/picard.jar`
`MarkDuplicates`
`INPUT=harismendy_data/class/P21.sorted`
`.bam OUTPUT=P21.sorted.rmdup.bam`
`METRICS_FILE=myrmdupMetrics.txt`
`REMOVE_DUPLICATES=true`
`ASSUME_SORTED=true`

Stats and Slices

Calculate flag statistics

```
samtools flagstat  
harismendy_data/class/P21.sorted.bam > P21.flagstat.txt
```

How many “gapped reads” ?

```
samtools view harismendy_data/class/P21.sorted.bam |  
awk '$6~/[ID]/' | wc -l
```

Subset the reads from an interval

```
samtools view -bh -L harismendy_data/class/TSG.bed  
harismendy_data/P21.sorted.bam > P21.sorted.TSG.bam
```

or

```
samtools view -bh harismendy_data/P21.sorted.bam  
chr1:120454175-120612317 > P21.NOTCH2.bam
```

Visualize Alignments in IGV

- Transfer NOTCH2 BAM to your laptop (WinSCP, Fugu, cyberDuck, scp)
- Start IGV
- Use human hg19 genome reference

BEDTOOLS coverage

Read the doc

<http://bedtools.readthedocs.org/en/latest/>

Read examples

<https://github.com/arg5x/bedtools-protocols/blob/master/bedtools.md>

Calculate coverage depth over CGC genes (chr1 only)

```
bedtools coverage -hist -abam  
harismendy_data/class/AA2253B_groupRealigned.chr1.bam -b  
harismendy_data/class/CGC.exons.chr1.bed >  
AA2253B.CGC.hist.cov.txt
```

Fraction of CGC bp covered at >30x ?

```
grep '^all' AA2253B.CGC.hist.cov.txt | awk '$2>30' | awk  
{sum=sum+$5} END {print sum}'
```

Calculate HS Metrics

- ```
java -jar /opt/biotools/picard/picard.jar
CalculateHsMetrics
BAIT_INTERVALS=harismendy_data/resources/humanV4-
baits_hg19_lite.interval_list
TARGET_INTERVALS=harismendy_data/resources/resource/huma
nV4-targets_hg19_lite.interval_list
INPUT=harismendy_data/class/P21.sorted.bam
OUTPUT=P21.HsMetrics.txt
```